

Catalog No	Product Name	Detection Range	Research Area
RE1060MF	MF-Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	10.94-700pg/mL	Cancer, Metabolism
RE1026MF	MF-Mouse GRO α /CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	15.63-1000pg/mL	Immunology, Neuroscience
RE1059MF	MF-Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63-1000pg/mL	Cancer, Immunology
RE3186HF	MF-Human IL-6(Interleukin 6) ELISA Kit	0.78-50pg/mL	Immunology,Stem cells
RE2839HF	MF-Human IFN- α (Interferon Alpha) ELISA Kit	3.90-250pg/mL	Cancer, Immunology
RE2029H	Human VF(Visfatin) ELISA Kit	0.31-20 ng/mL	metabolic disease
RE1094H	Human TIMP1(Tissue Inhibitors Of Metalloproteinase 1) ELISA Kit	0.16-10 ng/mL	hypercholesterolemia disease
RE2088H	Human LZM(Lysozyme) ELISA Kit	1.25-80 ng/mL	Nephropathy
RE2854H	Human CRP(C Reactive Protein) ELISA Kit	125.00-8000 pg/mL	inflammation、 cardiac disease
RE2846M	Mouse FABP1(FattyAcid Binding Protein 1, Liver) ELISA Kit	0.16-10ng/mL	metabolic disease、 cancer
RE3189M	Mouse POSTN/OSF-2(Periostin) ELISA Kit	0.16-10ng/mL	Cell & Developmental Biology
RE1737M	Mouse S100B(S100Calcium Binding Protein B) ELISA Kit	78.13-5000pg/mL	nerve system disease
RE2827M	Mouse LBP(Lipopolysaccharide Binding Protein) ELISA Kit	3.13-200ng/mL	inflammation
RE1979H	Human FBLN1(Fibulin 1) ELISA Kit	0.63-40 ng/mL	signal pathway、 Cell
RE1147M	Mouse FABP4(FattyAcid Binding Protein 4, Adipocyte) ELISA Kit	0.78-50ng/mL	inflammation
RE1730M	Mouse S100A9(S100Calcium Binding Protein A9) ELISA Kit	0.63-40ng/mL	Immunology、 neurology
RE3205M	Mouse PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit	62.5-4000pg/mL	hypercholesterolemia disease
RE2804H	Human APOB48(Apolipoprotein B48) ELISA Kit	1.56-100 ng/mL	hypercholesterolemia disease
RE2985M	Mouse SELE(E-selectin) ELISA Kit	0.32-20ng/mL	metabolic disease

TEL: +86 027 65316809
 E-mail: service@reedbiotech.com
 Website: www.reedbiotech.com
 Address: Building B3-3, Biolake, No. 666 Gaoxin Avenue, East Lake High-tech Zone, Wuhan, PRC


Reed Biotech Ltd



ELISA KIT
ANTIBODY
CREATIVE CENTER

High quality & High sensitivity

Micro samples, Fast speed tests



Universal series
Sample volume: 50μL
Suitable for all types of biological liquid samples
Only 90min test time



TECHNOLOGY ADVANTAGE

Core Reagents: Building consistent, reproducible assays hinges on high-quality raw materials. At the center of all Reed Biotech immunoassays are our proteins, antibodies, and specialized diluents.

Expertise: Our deep experience in developing and manufacturing immunoassays allows us to produce the most consistent and reliable assays capable of delivering superior accuracy and precision.

Accuracy: Assay interference common in complex sample matrixes can be overcome with careful choice of reagents and specialized diluents, helping to ensure you get the most accurate and reliable results.

Quality Control: Rigorous quality control of every raw material and every assay kit allows us to deliver assays that remain consistent and highly reproducible over time.

Custom Services: Build the assay that best fits for your experimental needs using our customized immunoassay solutions. Services include species testing, qualification, modification of an existing assay, and new immunoassay development.

Fast, stable, saving, shorter experiment time, 4°C storage temperature, 25μL loading volume breaking the Mouse limit (the same as some of the indicators in Human & Rat)



2-8°C
Bulk storage



25μL
Microloading



Micro & Fast
Roll into one

Cat. No	Product Name	Sensitivity	Detection Range	Research Area
RE10185	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	0.94 ng/mL	1.56-100 ng/mL	Immunology, Cancer
RE10165	MDA(Malondialdehyde) ELISA Kit	9.15 ng/mL	31.25-2000 ng/mL	Neuroscience
RE10133	ALD(Aldosterone) ELISA Kit	9.41 pg/mL	31.25-2000 pg/mL	Signal transduction
RE10132	NE(Noradrenaline) ELISA Kit	26.1 pg/mL	78.13-5000 pg/mL	Neuroscience, Cardiovascular
RE10092	cAMP(Cyclic Adenosine Monophosphate) ELISA Kit	92.5 pg/mL	312.5-20000 pg/mL	Signal transduction
RE10118	PG(Progesterone) ELISA Kit	0.47 ng/mL	1.57-100 ng/mL	Signal Transduction
RE10082	cGMP(Cyclic Guanosine Monophosphate) ELISA Kit	45.2 pg/mL	156.25-10000 pg/mL	Signal transduction
RE10131	PGE2(Prostaglandin E2) ELISA Kit	0.97 pg/mL	3.13-200 pg/mL	Signal Transduction
RE10135	DA(Dopamine) ELISA Kit	4.71 pg/mL	15.63-1000 pg/mL	Neuroscience
RE10166	5-HT(5-Hydroxytryptamine) ELISA Kit	3.97 ng/mL	14.07-900 ng/mL	Neuroscience
RE10134	EPI(Epinephrine) ELISA Kit	9.74 pg/mL	31.25-2000 pg/mL	Cell Biology, Metabolism,Cancer
RE10175	Testo(Testosterone) ELISA Kit	17.7 pg/mL	62.5-4000 pg/mL	Epigenetics and Nuclear Signal,
RE10008	GnRH(GonadotropinReleasing Hormone) ELISA Kit	9.38 pg/mL	15.63-1000 pg/mL	Developmental biology,
RE10013	TGF-β1(TransformingGrowth Factor Beta 1) ELISA Kit	0.10 ng/mL	0.16-10 ng/mL	Cancer, Cardiovascular,
RE10142	HA(Histamine) ELISA Kit	0.51 ng/mL	1.57-100 ng/mL	Immunology
RE10172	gABA(Gamma-Aminobutyric Acid) ELISA Kit	8.73 pg/mL	31.25-2000 pg/mL	Neuroscience
RE3203M	Mouse SELP(P-Selectin) ELISA Kit	0.1 ng/mL	0.16-10 ng/mL	Cancer, Cardiovascular, Immunology, Stem cells
RE1015M	Mouse MPlF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Immunology
RE1098M	Mouse C5a(Complement Component 5a) ELISA Kit	37.5 pg/mL	62.5-4000pg/mL	Immunology, Signal transduction
RE1765M	Mouse OPG(Osteoprotegerin) ELISA Kit	37.5 pg/mL	62.5-4000pg/mL	Cancer, Cardiovascular
RE1059M	Mouse IFN-γ(Interferon Gamma) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Cancer, Immunology, Stem cells
RE2867M	Mouse IL-17A(Interleukin 17A) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Cancer, Immunology
RE3187M	Mouse IL-10(Interleukin 10) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Cancer, Immunology
RE3186M	Mouse IL-6(Interleukin 6) ELISA Kit	18.75 pg/mL	31.25-2000pg/mL	Cancer, Metabolism, Immunology, Neuroscience
RE1055M	Mouse IL-4(Interleukin 4) ELISA Kit	18.75 pg/mL	31.25-2000pg/mL	Cancer, Immunology
RE2866M	Mouse IL-2(Interleukin 2) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Cancer, Immunology
RE1074M	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	4.69 pg/mL	7.82-500pg/mL	Cancer, Cardiovascular, Metabolism, Immunology
RE2862M	Mouse GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	9.38 pg/mL	15.63-1000pg/mL	Cancer, Immunology
RE3202M	Mouse RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	18.75 pg/mL	31.25-2000pg/mL	Cancer, Cardiovascular, Immunology
RE2831M	Mouse ICAM-1/CD54(Intercellular Adhesion Molecule 1) ELISA Kit	0.82 ng/mL	1.37-1000 ng/mL	Cancer, Cardiovascular, Immunology, Stem cells



Combination of stable, micro sample detection, shorten the experimental steps. It is one of the flagship products that Reed Biotech focuses on R&D and promotion. In addition to being used for conventional Mouse, experimental species such as Human and Rat are being developed simultaneously. For details, log in to <http://www.reedbiotech.com/list/109.html>